



# ADULT INTERMITTENT URINARY CATHETER DECISION AID

**Aseptic Non-Touch Technique (ANTT) must be followed**

A sterile, single use intermittent catheter required for each urinary catheterization as per Health Canada licence  
Nurses should adhere to organizational policies and procedures and nursing regulatory body

## Female Anatomy

Typical sizes:  
10 (black), 12 (white), 14 (green)

## Size

The smallest French (FR) size that enables insertion and adequate drainage should be selected. Funnel end has colour-coding to allow for ease of recognition

## Male Anatomy

Typical sizes:  
14 (green), 16 (orange), 18 (red)

## Female Anatomy

Typical length:  
7-22 cm (may use male length for ease of drainage into a receptacle)

## Length

Male anatomy requires a longer length intermittent catheter due to longer urethra

## Male Anatomy

Typical length:  
33-45 cm (do not use female length)

## Hydrophilic (preferred)

- PU, PVC, silicone, and TPE
- activation of lubricated surface can be ready-to-use or require a package to be broken
- no added lubrication needed
- reduces friction and trauma
- may reduce the risk of UTI

## Material

Latex not recommended due to latex allergies and increasing sensitivities

## Prelubricated

- PVC, silicone, and TPE
- sterile water-soluble lubricant is already added to an uncoated intermittent catheter

## Uncoated

- PVC, silicone, and TPE
- sterile single use packet of water-soluble lubricant must be added to surface

## Smooth eyelets (preferred)

- fire polished or ultrasonically smoothed
- minimize urethral abrasion and erosion

## Eyelets

## Nonsmoothed eyelets

- cold punched eyelets
- can cause urethral abrasion and erosion
- may cause insertion pain and bleeding

## Standard/Straight/Nelaton

- common
- straight tip
- used when no anticipated challenges with catheterization

## Coudé/Tiemann/Curved

- curved firm tip
- insert using guide indicator
- often used for an enlarged prostate or urethral stricture

## Tips

## Rounded/Ergothan

- flexible bead tip allowing for easier catheterization
- often used for enlarged prostate, urethral strictures, or spastic pelvic floor

## Olive

- rounded, larger bulb tip
- often used for urethral strictures or females having difficulty locating the urethral meatus

## Standard

- basic intermittent catheter

## Compact

- smaller package
- may be easily carried in pocket or purse
- discreet

## System

## Closed System

- intermittent catheter and drainage bag is an all-in-one system
- useful when toilet access is limited
- urine volume indicator

## Protective Sleeve or Grip

- touch-less / no-touch
- decreases insertion contamination caused by handling
- useful when user needs to touch their catheter or when learning

## Financial

- consider public and private insurance, where available provincial/territorial/federal reimbursement programs
- consider patient support programs through manufacturers to navigate coverage

## Additional Considerations

## Environmental

- contact manufacturers for more details
- consult your local waste management policy for disposal of intermittent catheter and packaging

## Clinical

### Poor Hand Dexterity

- consider touchless with or without introducer tip
- consider use of insertion aids/devices

### Positioning

- females may find shorter catheters are easier to grasp and insert while seated on toilet
- females can use longer catheters to empty their bladder into a toilet or container while in a wheelchair

### Urine Volume Measurement Required

- consider use of closed system (with attached bag containing a volume scale)
- may use a urine collection container with a volume scale to drain and measure output

## Self-Management Considerations

### Discretion

- consider discrete/compact type
- consider use of closed system with attached catheter and urine bag attached if traveling or unable to access a toilet or appropriate facility
- consider discreet packaging
- consider discreet disposal features

## Social

### Financial Limitations

- inform patient of catheter cost ranges to self-determine choice
- consider individual's financial limitations when selecting intermittent catheters after insurance and reimbursement options have been confirmed
- consider patient support programs provided by manufacturers and distributors

### Education

- evaluate patient's health literacy and readiness for self-care with each clinical encounter

## Insertion

- practice ANTT
- *always* do pericare prior to any catheterization
- when doing pericare examine urethra for any aberrations that may affect type of tip selection
- review chart prior to selection to note any conditions that may affect ease of insertion (e.g., prolapse, obstruction, stricture)
- select the most appropriate catheter for the individual and document any special considerations to ensure continuity
- consider using local topical anesthetic application prior to catheterization if the patient may be fearful or has sensitivities

## Removal

- *slowly* withdraw catheter 2.5 cm (1") at a time, rotating and repositioning each time you withdraw to ensure complete bladder emptying

## Troubleshooting

- with no expected urine return, the catheter may be curled inside urethra, so use a firmer catheter, Coudé tip, or rounded/Ergothan tip
- bladder scan post catheterization to check complete emptying, if available

**Notes:** Abbreviations PU = polyurethane, PVC = polyvinyl chloride, TPE = thermoplastic elastomer, UTI = urinary tract infection  
 Nurses must practice within nursing regulatory body scope of practice and in accordance with health care organizational policies and procedures. Always refer to manufacturer's instructions for use.  
 Developed by a nursing task force through Nurses Specialized in Wound, Ostomy and Continence Canada (NSWOCC). Sponsored by an unrestricted educational grant from Coloplast Canada. Trademarks acknowledged.